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
AASM, Clinical Guideline for the Evaluation and Management of Chronic Insomnia in Adults; CBT-I, cognitive behavior therapy for insomnia; ESRS, European Sleep Research Society

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Study design This retrospective cohort study was conducted using longitudinal data collected via a mobile app named Good Sleep 365 between November 14, 2018, and February 28, 2022. Three therapeutic modes (ie, dCBT-I, medication, and their combination) were compared at month 1, month 3, and month 6 (primary). Inverse probability of treatment weighting (IPTW) using propensity scores was applied to enable homogeneous comparisons between the 3 groups. EXPOSURES Treatment with dCBT-I, medication therapy, or combination therapy according to prescriptions. The Pittsburgh Sleep Quality Index (PSQI) score and its essential subitems were used as the primary outcomes. A total of 4052 patients (mean [SD] age, 44.29 [12.01] years; 3028 [74.7%] female) participants were selected for dCBT-I (n = 418), medication (n = 862), and their combination (n = 2772).


CI, confidence interval; dCBT-I, digital cognitive behavior therapy for insomnia; PSQI, Pittsburgh Sleep Quality Index

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
Self-report Sleep metrics before and after CBT-I

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Post-Treatment Medication Use by Medication Category

Study design To evaluate the effectiveness, feasibility, and acceptability of a telehealth group-based CBT-I program that incorporates gradual dose reduction (GDR) protocols for sleep-related medication and substance use. We conducted a retrospective cohort study of 204 adults with DSM-5-diagnosed insomnia disorder who participated in a 7-week CBT-I group medical visit program at an academic integrative health clinic. The program included one intake session, five weekly group sessions, and one follow-up session. Individualized GDR protocols targeted a broad range of substances. We assessed pre-post changes in insomnia symptoms, consensus sleep diary (CSD) metrics, medication/substance use, and treatment acceptability.

CBD, cannabidiol; CBN, cannabinol; CBT-I, cognitive behavior therapy for insomnia; GDR, gradual dose reduction; THC, tetrahydrocannabinol

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
Primary endpoint: Sleep efficiency

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
Change from Baseline of SE at week 7 (FAS)

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Change from Baseline of DBAS-16 at week 7 (FAS)

Study design 68 adults diagnosed with insomnia were enrolled in a double-blind, sham-controlled randomized decentralized clinical trial. All participants completed one on-site visit for written consent, after which all procedures—including intervention use, sleep diary recording, and questionnaire assessments—were conducted remotely over 6-7 weeks. Sleep efficiency (primary outcome) and other sleep-related measures were collected through the app. Group differences were analyzed using least-squares estimates.

DBAS-16, Dysfunctional Beliefs and Attitudes about Sleep - 16 items; SE, Sleep Efficiency; ISI, Insomnia severity index; LSM, Least Squares Mean; se, standard error; FAS, Full Analysis Set

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References. 1. 대한신경정신의학회. 2019 불면증 임상진료지침. 불면증의 진단과 치료 2. 식품의약품안전처. 의료기기 제조 허가증. 제허 23-511호 3. McNamara S, et al. Chronic Insomnia. Available at: <https://www.ncbi.nlm.nih.gov/books/NBK526136/> (accessed on 17 Apr. 2026) 4. Lee E. Cognitive behavior therapy for insomnia. J Korean Med Assoc. 2020;63(8):477-483 5. Riemann D, et al. The European Insomnia Guideline: An update on the diagnosis and treatment of insomnia 2023. J Sleep Res. 2023;32(6):e14035 6. Sateia MJ, et al. Clinical Practice Guideline for the Pharmacologic Treatment of Chronic Insomnia in Adults: An American Academy of Sleep Medicine Clinical Practice Guideline. J Clin Sleep Med. 2017;13(2):307-349 7. Lu M, et al. Comparative Effectiveness of Digital Cognitive Behavioral Therapy vs Medication Therapy Among Patients With Insomnia. JAMA Netw Open. 2023;6(4):e237597 8. Gutierrez G, et al. Outcomes from a combined cognitive behavioral therapy for insomnia (CBT-I) and sleep-related medication and substance use reduction treatment. Sleep Med. 2026;138:108696 9. European Commission. EUDAMED-European Database on Medical Devices > UDI/Devices > SleepQ. Available at: <https://ec.europa.eu/tools/eudamed/#/screen/search-device?sm=KR-MF-000040293&deviceStatusCode=refdata.device-model-status.on-the-market&paging=%7B%22pageIndex%22%3A0%22pageSize%22%3A25%7D&sorting=%7B%22sortField%22%3AprimaryDi%22%22sortDirection%22%3Adesc%22%7D&submitted=true> (accessed on 09 Jul. 2026) 10. KSR ISO/IEC 27001 certificate (인증번호: ISMC-0032) 11. Park KM, et al. Validating the Efficacy of a Mobile Digital Therapeutic for Insomnia (WELT-I): Randomized Controlled Decentralized Clinical Trial. J Med Internet Res. 2025;27:e70722 12. 보건복지부 고시 제2025-6호 13. Thakral M, et al. ISI-3: Evaluation of a brief screening tool for insomnia. Sleep Med. 2021;82:104-109 14. Kim M, et al. Validation of the Patient Health Questionnaire-9 and Patient Health Questionnaire-2 in the General Korean Population. Psychiatry Investig. 2023;20(9):853-860 15. Cho A, et al. The Effects of Generalized Anxiety Disorder 7-Item Scale on Health Behavior among Korean Adults: The 8th Korea National Health and Nutrition Examination Survey (2021). Korean J Fam Pract. 2024;14(2):75-81 16. Scharf MT. Reliability and Efficacy of the Epworth Sleepiness Scale: Is There Still a Place for It? Nature and Science of Sleep 2022;14 2151-2156 17. Morin CM, et al. Dysfunctional Beliefs and Attitudes about Sleep (DBAS): Validation of a Brief Version (DBAS-16). Sleep. 2007;30(11):1547-1554 18. Ong AA, Gillespie MB. Overview of smartphone applications for sleep analysis. World J Otorhinolaryngol Head Neck Surg. 2016;2(1):45-49 19. Basiouny DE, Habib HI. Cognitive Behavioral Therapy for Insomnia (CBT-I): A Comprehensive Review. Indian Journal of Sleep Medicine. 2023;18(4):62-65 20. Hwang JW, et al. Systematic review and meta-analysis on fully automated digital cognitive behavioral therapy for insomnia NPJ Digit Med. 2025;8(1):157 21. Mitchell MD, et al. Comparative effectiveness of cognitive behavioral therapy for insomnia: a systematic review. BMC Fam Pract. 2012;13:40

Selected Product Information [의료기기]

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